pNAPEP-1025

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Version 2022/11

pNAPEP-1025

Activated factor X substrate
(identical to CBS 3139™)

REF 61031025

RUO

Research Use Only

CRYOPEP

83 rue Yves Montand
F-34080 MONTPELLIER
Tél.: +33(0)4 67 10 71 20
Fax : +33(0)4 67 10 71 21
contact@cryopep.com
www.cryopep.com

COMPOSITION

Vial : 25 mg pNAPEP-1025, 0.0414 mmol
Bulking agent : (D)-mannitol (80 mg/vial)

The accurate quantity is indicated on analytical data sheet.

APPLICATIONS

This chromogenic substrate could be used by all methods using activated factor X (FXa).

REAGENT PREPARATION

Reconstitute the vial with 20 mL of sterile distilled water to obtain a concentration around 2 mmol/L.

Let stand at room temperature (18-25°C) for 30 minutes.

Shake before use.

CHEMICAL STRUCTURE

C₂₁H₃₄N₈O₇S, AcOH

Formula : CH₃SO₂-D-Leu-Gly-Arg-pNA,AcOH
(Methanesulfonyl-D-leucyl-glycyl-L-arginine-paranitroaniline acetate)

Molecular weight : 602.66 g/mol

PHYSICO-CHEMICAL DATA

Solubility : >10 mmol/L (H₂O)

ε_{316 nm} : 1,27 · 10⁴ mole⁻¹ · L · cm⁻¹

ε_{405 nm} : 1 · 10⁴ mole⁻¹ · L · cm⁻¹

STABILITY

Lyophilized powder :

Stable until the expiry date, if stored at 2-8°C in the dark. Avoid light exposure.

Keep in a dry place, hygroscopic product. Shipment of product does not require cooling during the time of transportation.


Aqueous solution (1 - 5 mmol/L) :

Stable for 6 months at 2-8°C.

PRECAUTIONS

- The product should be sealed in its original vial, following each use to ensure its stability.
- Minimize the evaporation of product during use.
- Contamination by microorganism may cause hydrolysis.
- Yellowing of the substrate indicates contamination, discard and use a new vial.

PRINCIPLE

The rate of pNA, released by enzymatic hydrolysis, is quantified by spectrophotometry at 405 nm and is proportional to the enzyme activity

ENZYMATIC DATA

Bovine Factor Xa : K_m = 2,2 · 10⁻⁴ mol/L, k_{cat} = 300 s⁻¹

(These data are determined in Tris HCl 50 mmol/L ; NaCl 110 mmol/L ; buffer pH 8.3 at 37°C)

Version 2022/11