pNAPEP-1751

Substrat de la plasmine
(identique au S-2251TM de Chromogenix)

REF 6101-1751

RUO

Usage Recherche

CRYOPEP

83 rue Yves Montand
F-34080 MONTPELLIER
Tél.: +33(0)4 67 10 71 20
Fax : +33(0)4 67 10 71 21
contact@cryopep.com
www.cryopep.fr

COMPOSITION

Flacon de 25 mg de pNAPEP-1751, soit 0.0453 mmol

Excipient : (D)-mannitol (80 mg/flacon)

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pNAPEP-1751

Plasmin substrate
(identical to S-2251TM from Chromogenix)

REF 6101-1751

RUO

Research Use Only

CRYOPEP

83 rue Yves Montand
F-34080 MONTPELLIER
Tél.: +33(0)4 67 10 71 20
Fax : +33(0)4 67 10 71 21
contact@cryopep.com
www.cryopep.com

COMPOSITION

Vial : 25 mg pNAPEP-1751, 0.0453 mmol

Bulking agent : (D)-mannitol (80 mg/vial)

The accurate quantity is indicated on analytical data sheet.

APPLICATION

The substrate has been used for plasmin and. This chromogenic substrate could be used for the determination in plasma of :

- Antiplasmin in plasma
- Plasmin in plasma
- streptokinase-activated plasminogen

REAGENT PREPARATION

Reconstitute the vial with 20 mL of sterile distilled water to obtain a concentration up to 2.26 mmol/L.

Let stand at room temperature (18-25°C) for 30 minutes.

Shake before use.

CHEMICAL STRUCTURE

C₂₃H₃₈N₆O₅, HCl

Formula : H-D-Val-Leu-Lys-pNA, 2 HCl
(H-D-Valyl-L-Leucyl-L-Lysine-paranitroaniline dihydrochloride)

Molecular weight : 551,49 g/mol

PHYSICO-CHEMICAL DATA

Solubility : > 40 mmol/L (H₂O)

$\epsilon_{316\text{ nm}}$: $1,27 \cdot 10^4 \text{ mole}^{-1} \cdot \text{L} \cdot \text{cm}^{-1}$

STABILITY

Lyophilized powder :

Stable until the expiry date, if stored at 2-8°C.

Avoid light exposure.

Keep in a dry place, hygroscopic product.


Aqueous solution (3 - 4 mmol/L) :

Stable for 6 months at 2-8°C.

PRECAUTIONS

- The product should be sealed in its original vial, following each use to ensure its stability.
- Minimize the evaporation of product during use.
- Contamination by microorganism may cause hydrolysis.
- Yellowing of the substrate indicates contamination, discard and use a new vial.

PRINCIPLE

The rate of pNA, released by enzymatic hydrolysis, is quantified by spectrophotometry at 405 nm and is proportional to the enzyme activity

ENZYMATIC DATA

The following kinetic constants were estimated at 37°C in 2,5 mL 0.05 mol/L Tris buffer pH 7.4, I 0.5.

Enzyme

Human Plasmin: $K_m = 3 \cdot 10^{-4} \text{ mol/L}$;
 $V = 0,5 \times 10^{-6} \text{ mol/min CU}$

Plasminogen SK: $K_m = 2 \cdot 10^{-4} \text{ mol/L}$;
 $V = 1 \times 10^{-6} \text{ mol/min CU}$

SELECTIVITY

- This substrate is insensitive to kallikrein (glandular and plasma) and urokinase.