

DEFICIENT PLASMAS

Immunodepleted deficient plasmas

Fresh Frozen Plasmas

Prekallikrein Immunodepleted Deficient Human Plasma



Associated products

a2-Antiplasmin Immunodepleted Deficient Human Plasma

Antithrombin deficient plasma immuno depleted

Antithrombin Immunodepleted Deficient Human Plasma

Informations

Prekallikrein is a glycoprotein, zymogen of serine protease. Non-covalently complexed with high molecular weight kininogen.

Prekallikrein participates in the activation of coagulation, fibrinolysis, generation of kinins and inflammatory phenomena.

It is activated into kallikrein by FXIIa.

Reference	Presentation	Format
6-FDPPK	Bottle	1 x 100 mL
6-FDPPK-10	Kit	10 x 1.0 mL

Immunodepleted deficient plasma for the determination of prekallikrein.

Citrated normal human plasma depleted in prekallikrein by antibodies specific to prekallikrein grafted on agarose gels. Contains 20 mM Hepes buffer.

Components

- 10 cryotubes x 1 mL or 100 mL vial of frozen plasma.

Advantages

- No bovine additives
- No reconstitution error
- No plasma alteration linked to freeze-drying
- Ready to use after defrosting (4 min at 37° C)
- Cryotubes ready to use after thawing (4 min at 37° C).

Characteristics

Packaging in plastic cryotubes or in bottles of at least 100 mL.

Frozen, immuno-depleted plasmas are certified to have less than 1% for the deficient factor considered, both for the antigenic and functional assay in haemostasis.

This box is intended for research use