

HUMAN PLASMAS

Plasmas and interferences

Plasma with ultra high triglycerides level



Associated products

Hemolysed Plasma

Icteric plasma with high bilirubin

Lipemic plasma

Informations

Triglycerides are one of the lipid compounds essential to the body. They allow the storage of fatty acids in adipose tissue and come either from the diet or be synthesized by the body. An excess of triglycerides (hypertriglyceridemia) constitutes a cardiovascular risk factor.

Advantages

Minimize test time. Ready to use.

Characteristics

Special plasmas are derived from patients with a congenital deficiency, severe or moderate, or presenting a particular profile. No buffer or preservatives are added. Quickly frozen at -80°C , the plasma maintains perfectly intact the matrix. All plasmas are stable when stored at -40°C to -80°C . We carefully pack with dry ice during shipment. No additive or preservative. Expiry date > 1 year. Plastic vials.

