

POLYCLONAL ANTIBODIES

Anti-heparin

Sheep polyclonal antibody anti-Human heparin coFactor II



Informations

The second heparin cofactor is a serine protease inhibitor. It inhibits thrombin, chymotrypsin and other enzymes of the same group. Its rate of inhibition is amplified in the presence of heparin.

Reference	Presentation	Format
9-PAHCII-S	Vial	1 mg

Antigen: human heparin cofactor II

Application : Immunoblotting, ELISA

Host : Sheep

Immunogen: Purified human heparin cofactor II

Advantages

Custom needs by supplying you antibodies conjugated with biotin, HRP, FITC or other conjugates. Special formulations are available upon request. Discount according to quantities

Characteristics

The vast majority of antibodies is pure (without additives) with > 95 % purity SDS-PAGE. Stock antibodies are supplied in 50 % glycerol/water (v/v) for ease of storage and use. Both small, laboratory scale and bulk, production scale quantities are available. Expiration date of one year from delivery.

