

# POLYCLONAL ANTIBODIES

## Anti-plasminogen

### Sheep polyclonal antibody anti-Human plasminogen



#### Associated products

Sheep pAb anti-mouse plasminogen

#### Informations

Plasminogen is the zymogen of plasmin, a key enzyme in the fibrinolysis system. Plasminogen is synthesized mainly by the liver but also the eosinophils, the kidney and the cornea. It exists in 2 molecular forms: glu-plasminogen (native form) and lys-plasminogen (more active form). The main pathways for activating plasminogen to plasmin involve t-PA and u-PA.

Reference	Presentation	Format
9-PAHPG-S	Vial	1 mg

#### Antigen: human plasminogen

Application : Immunoblotting, ELISA

Host : Sheep

Immunogen : Purified human plasminogen

#### Advantages

Custom needs by supplying you antibodies conjugated with biotin, HRP, FITC or other conjugates. Special formulations are available upon request. Discount according to quantities

#### Characteristics

The vast majority of antibodies is pure (without additives) with > 95 % purity SDS-PAGE. Stock antibodies are supplied in 50 % glycerol/water (v/v) for ease of storage and use. Both small, laboratory scale and bulk, production scale quantities are available. Expiration date of one year from delivery.

