

ANTICOAGULANT MONITORING

ANTI-Xa

ORGARAN®

ANTICOAGULANT CONTROLS

Lyophilized plasmas

TECHNOVIEW® Orgaran® Cont High



Associated products

TECHNOVIEW® Orgaran® Cal Set

TECHNOVIEW® Orgaran® Cont Low

Informations

Danaparoid sodium (Orgaran®) is a polysaccharide anticoagulant, used as an alternative treatment to unfractionated heparins and low molecular weight heparins when the latter are contraindicated.

Reference	Presentation	Format
4-5090114	Vial	5 x 1.0 mL

High control plasmas for the Orgaran® assay.

TECHNOVIEW® Orgaran® Control High quality control plasmas are titrated to approximately 1.0 IU / mL, optimized for anti-FXa methods.

Components

- 5 vials x 1 mL lyophilized plasma

Characteristics

TECHNOVIEW® Orgaran® plasmas are prepared from citrated plasmas supplemented with different concentrations of sodium danaparoid. They make it possible to validate the sodium danaparoid calibration curve.

