

# ANTICOAGULANT MONITORING

DOAC

DABIGATRAN

ANTICOAGULANT CONTROLS

Chronometric assay

## TECHNOVIEW® Dabigatran Cont Low



### Associated products

TECHNOVIEW® Dabigatran Cal Set

TECHNOVIEW® Dabigatran Cont High

### Informations

Dabigatran is a potent, direct, competitive and reversible thrombin inhibitor. It inhibits free and fibrin-bound thrombin and thrombin-induced platelet aggregation.

Reference	Presentation	Format
4-5090214	Vial	5 x 1.0 mL

### Low control plasmas for dabigatran assay

TECHNOVIEW® Dabigatran Low Control quality control plasmas are prepared from citrated supplemented plasmas titrated to approximately 100 ng / mL, optimized for anti-FIIa methods.

### Components

- 5 vials X 1 mL lyophilized plasma

### Characteristics

TECHNOVIEW® Dabigatran plasmas are prepared from citrated plasmas supplemented with different concentrations of dabigatran.

