

LYOPHILIZED CALIBRATORS AND CONTROLS

SCREENING TEST CONTROLS

MULTIPARAMETRIC CONTROLS

Lyophilized plasmas

TECHNOCLOT® Control A



Associated products

DAPTIN® TC
Fibrinogen Reagent
Siron LIS (aPTT liquid)
Siron LS (aPTT liquid)
TECHNOPLASTIN® HIS

Reference	Presentation	Format
4-5021070	Vial	10 x 1.0 mL
4-5021075	Vial	50 x 1.0 mL

Control plasma between normal and abnormal activity in hemostasis.

This quality control is titrated for routine hemostasis tests (PT, aPTT, TT and fibrinogen) for abnormal values.

Informations

PT : Prothrombin level
aPTT : Activated partial thromboplastin time
TT : Thrombin time

Components

- 10 or 50 vials x 1 mL lyophilized plasma

Characteristics

TECHNOCLOT® Control A is an abnormal citrated human plasma in which the level of coagulation factors is reduced.
TECHNOCLOT® Control A contains stabilizers but no bactericidal adjuvants.

