

**CERTIFICAT D'ANALYSES**


**CRYOcheck™ Heparin Control**

Référence  
**CCH-10**

Lot  
....

Péremption  
....

**EXCLUSIVEMENT réservé au Diagnostic In Vitro**



Les résultats ont été obtenus dans les conditions opératoires suivantes :

- Instrument : STA-R®
- Réactif TCA : HemosIL APTT-SP

Pour utiliser ce plasma en tant que contrôle, chaque laboratoire doit établir ses propres valeurs à partir d'au moins 20 essais en séries indépendantes (Cembrowski GS, Carey RN. *Laboratory quality management*. Chicago: ASCP Press; 1989).

Les résultats ci-dessous sont calculés à partir de données spécifiques obtenues à partir d'une combinaison de réactifs et d'instruments utilisés pour des tests de contrôle de qualité interne. Ils peuvent varier selon l'instrument, la méthode et la technique utilisés.

	TCA
Moyenne	92,3 s

Contrôle Qualité

**Fiche de sécurité**

Disponible sur simple demande auprès de notre service commercial



**Precision Biologic Inc.**  
140 Eileen Stubbs Avenue  
Darmouth B3B 0A9, Nova Scotia, Canada



**Emergo Europe**  
Prinsessegracht 20  
2514 AP The Hague, Pays-Bas



**CRYOPEP**  
83, rue Yves Montand  
34080 MONTPELLIER

**CERTIFICATE OF ANALYSIS**

**CRYOcheck™ Heparin Control**

Catalog number  
**CCH-10**

Lot  
....

Expiry  
....

**ONLY for In Vitro Diagnostic**



Test system information :

- Instrument :
- APTT reagent :

STAR®  
HemosIL APTT-SP

For controls to work effectively, control ranges must be established in each laboratory. Optimally, control data should be collected from at least 20 runs on different days before control ranges are calculated (*Cembrowski GS, Carey RN. Laboratory quality management. Chicago: ASCP Press; 1989*).

The results presented below are calculated from data specific to the reagent-instrument combination used for internal quality control testing at Precision Biologic.

Results may vary depending on instrument, method, and technique used.

	APTT
Mean	92,3 s


*Quality Control*

**Material safety data sheet**

Available on request from our commercial department



**Precision Biologic Inc.**  
140 Eileen Stubbs Avenue  
Darmouth B3B 0A9, Nova Scotia, Canada



**Emergo Europe**  
Prinsessegracht 20  
2514 AP The Hague, Pays-Bas



**CRYOPEP**  
83, rue Yves Montand  
34080 MONTPELLIER