

CERTIFICAT D'ANALYSES

pNAPEP 0968
Lyophilisat

Référence
61030968

Lot
09681014

Péremption
2024/09

Usage RECHERCHE Exclusivement

RUO

2°C / 8°C

Conditions de stockage : Entre +2/+8°C à l'abri de la lumière et de l'humidité
Type : Substrat du facteur IXa
Séquence : H-D-Leu-PHG-Arg-pNA, 2AcOH
Poids moléculaire : 660,8 g/mol

	<i>Spécification</i>	Résultat	Conformité
Aspect	<i>Poudre blanche à jaune pâle</i>	Poudre blanche	Conforme
pH	<i>5 < < 8</i>	6,67	Conforme
pNA libre	<i>< 0,5 %</i>	0,28 %	Conforme
Teneur en peptide	<i>≥ 25 mg</i>	27.34 mg	Conforme
Reconstitution	<i>2 mmol/L</i>	20.69 mL d'eau distillée	

Contrôle Qualité



CERTIFICATE OF ANALYSIS

pNAPEP 0968
Lyophilisate

Catalog number
61030968

Lot
09681014

Expiry
2024/09

RESEARCH USE ONLY

RUO

2°C / 8°C

Storage conditions : +2/+8°C avoid light exposure and moisture
Type : Activated Factor IX substrate
Identity : H-D-Leu-PHG-Arg-pNA, 2AcOH
Molecular weight : 660,8 g/mol

	<u>Specification</u>	Result	Conformity
Appearance	<i>White to pale yellow powder</i>	White powder	Conform
pH	<i>5 < < 8</i>	6,67	Conform
Free pNA	<i>< 0,5 %</i>	0,28 %	Conform
Quantity /vial	<i>≥ 25 mg</i>	27.34 mg	Conform
Reconstitution	<i>2 mmol/L</i>	20.69 mL of distilled water	

Quality Control

