

CERTIFICAT D'ANALYSES

CRYOcheck™ Low Fibrinogen Control

Référence
CCLF-10

Lot
....

Péremption
....

EXCLUSIVEMENT réservé au Diagnostic In Vitro



Les résultats ont été obtenus dans les conditions opératoires suivantes :

- Instrument : Diagnostica Stago STAR®
- Réactif : Diagnostica Stago STA® Fibrinogen®

Pour utiliser ces contrôles de façon optimale, chaque laboratoire doit déterminer son intervalle de confiance à partir d'au moins 20 essais en séries indépendantes (*Cembrowski GS, Carey RN. Laboratory quality management. Chicago: ASCP Press; 1989*).

Les résultats ci-dessous sont calculés à partir de données spécifiques obtenues à partir d'une combinaison réactif / instrument utilisés pour des tests de contrôle de qualité interne. Ils peuvent varier selon l'instrument, la méthode et la technique utilisés.

Paramètre	Mesure	Méthode
Fibrinogène	0,93 g/L	Clauss

Contrôle Qualité

Fiche de sécurité

Disponible sur simple demande auprès de notre service commercial



Precision Biologic Inc.
140 Eileen Stubbs Avenue
Darmouth B3B 0A9, Nova Scotia, Canada



Emergo Europe
Prinsessegracht 20
2514 AP The Hague, Pays-Bas



CRYOPEP
83, rue Yves Montand
34080 MONTPELLIER

CERTIFICATE OF ANALYSIS

CRYOcheck™ Low Fibrinogen Control

Catalog number
CCLF-10

Lot
....

Expiry
....

ONLY for In Vitro Diagnostic



Test system information :

- Instrument : Diagnostica Stago STAR®
- Reagent : Diagnostica Stago STA® Fibrinogen®

For controls to work effectively, control ranges must be established in each laboratory. Optimally, control data should be collected from at least 20 runs on different days before control ranges are calculated (*Cembrowski GS, Carey RN. Laboratory quality management. Chicago: ASCP Press; 1989*).

The results presented below are calculated from data specific to the reagent/instrument combination used for internal quality control testing at Precision Biologic. Results may vary depending on instrument, method, and technique used.

Parameter	Value	Method
Fibrinogen	0,93 g/L	Clauss

Quality Control

Material safety data sheet

Available on request from our commercial department



Precision Biologic Inc.
140 Eileen Stubbs Avenue
Darmouth B3B 0A9, Nova Scotia, Canada



Emergo Europe
Prinsessegracht 20
2514 AP The Hague, Pays-Bas



CRYOPEP
83, rue Yves Montand
34080 MONTPELLIER