

HUMAN PLASMAS

High Factor plasmas

Plasma with high level of S protein: > 150 %



Associated products

Protein S human deficient plasma (acquired)

Protein S human deficient plasma (congenital)

Informations

Protein S is a vitamin K dependent protein. It is a physiological inhibitor of coagulation. It acts as a cofactor of activated protein C by promoting the inactivation of FVa and FVIIIa, prothrombin, of the prothrombinase complex, FX. A protein S deficiency can be either acquired (hepatocellular insufficiency, vitamin K deficiency, anti-protein S antibody, ...) or constitutional (heterozygous or homozygous deficiency) grouped into 2 types depending on whether the deficiency is quantitative (type I) or qualitative (type II).

Advantages

Minimize test time. Ready to use.

Characteristics

Special plasmas are derived from patients with a congenital deficiency, severe or moderate, or presenting a particular profile. No buffer or preservatives are added. Quickly frozen at -80° C, the plasma maintains perfectly intact the matrix. All plasmas are stable when stored at -40° C to -80° C. We carefully pack with dry ice during shipment. No additive or preservative. Expiry date > 1 year. Plastic vials.

