

HUMAN PLASMAS

Plasma with anti-phospholipid antibodies

Positive plasma to anti-prothrombin antibodies



Associated products

Positive plasma to other antiphospholipids antibodies

Informations

Anti-prothrombin antibodies are part of the anti-phospholipid antibody family. The anti-prothrombin antibodies would cause a transient deficit in prothrombin by binding to the prothrombin and by accelerating the elimination of the antigen-antibody complexes thus formed. About 20-30% of lupus have anti-prothrombin IgGs.

Reference	Presentation	Format
6-PPAII	Vial	1 x 1.0 mL

Advantages

Minimize test time. Ready to use.

Characteristics

Special plasmas are derived from patients with a congenital deficiency, severe or moderate, or presenting a particular profile. No buffer or preservatives are added. Quickly frozen at -80°C , the plasma maintains perfectly intact the matrix. All plasmas are stable when stored at -40°C to -80°C . We carefully pack with dry ice during shipment. No additive or preservative. Expiry date > 1 year. Plastic vials.

