

Plasma of insufficient kidney patient



Associated products

Hemolysed Plasma

Icteric plasma with high bilirubin

Lipemic plasma

Informations

Renal failure corresponds to a dysfunction of the kidneys which then no longer respond to their basic functions: production of hormones, maintenance of the acid-base balance, maintenance of the balance of minerals, filtration and purification of the blood and regulation of blood pressure by the production of hormones (renin and bradykinin).

Advantages

Minimize test time. Ready to use.

Characteristics

Special plasmas are derived from patients with a congenital deficiency, severe or moderate, or presenting a particular profile. No buffer or preservatives are added. Quickly frozen at -80°C , the plasma maintains perfectly intact the matrix. All plasmas are stable when stored at -40°C to -80°C . We carefully pack with dry ice during shipment. No additive or preservative. Expiry date > 1 year. Plastic vials.

Reference	Presentation	Format
6-PPIR	Vial	1 x 1.0 mL

