

HUMAN PLASMAS

Plasmas with anticoagulant drugs

Plasma with low molecular weight heparin (Innohep)



Associated products

Anticoagulant plasma – DTI – Argatroban – U/mL

Anticoagulant plasma – DTI – bivalirudin – U/mL

Plasma with direct thrombin inhibitor (Lepirudin)

Informations

Innohep is an anticoagulant from the low molecular weight heparin family. It prevents the formation or extension of clots in the blood vessels.

At low doses it is used in the preventive treatment of thromboembolic accidents in recent patients and at high doses, it is used in the treatment of venous thrombosis and pulmonary embolism.

Advantages

Minimize test time.
Ready to use.

Characteristics

Special plasmas are derived from patients with a congenital deficiency, severe or moderate, or presenting a particular profile.
No buffer or preservatives are added.
Quickly frozen at -80° C, the plasma maintains perfectly intact the matrix.
All plasmas are stable when stored at -40° C to -80° C.
We carefully pack with dry ice during shipment.
No additive or preservative.
Expiry date > 1 year.
Plastic vials.

