

ANTICOAGULANT MONITORING

DOAC

RIVAROXABAN

ANTICOAGULANT CONTROLS

Lyophilized plasmas



TECHNOVIEW® Rivaroxaban Cont High



Associated products

TECHNOVIEW® Rivaroxaban Cont Low

TECHNOVIEW® Rivaroxaban Cont Medium

Informations

Rivaroxaban is a highly selective direct inhibitor of FXa.

This inhibition interrupts the intrinsic pathway of the blood coagulation cascade, inhibiting the formation of thrombin and the development of thrombus.

It does not inhibit thrombin and has no effect on platelets.

Reference	Presentation	Format
4-5090174	Vial	5 x 1.0 mL

High control plasmas for the assay of rivaroxaban.

TECHNOVIEW® Rivaroxaban Control High quality control plasmas are titrated to approximately 300 ng / mL, optimized for anti-FXa methods.

Components

- 5 vials x 1 mL lyophilized plasma

Characteristics

TECHNOVIEW® Rivaroxaban (Xarelto®) plasmas are prepared from citrated plasmas supplemented with different concentrations of rivaroxaban (Xarelto®) and calibrated by HPLC / MS - MS method.

