

INSTRUMENTS

T-TAS®01

CONSUMABLES HD CHIPS - T-TAS® 01

ADDITIONAL REAGENTS

Anti-sedimentation reagent for HD Chip T-TAS®01



Associated products

HD Chip T-TAS® 01

Informations

Erythrocyte sedimentation inside the reservoir of the HD assay has been observed in some cases, and the prevention of erythrocyte sedimentation is expected to result in a more accurate and repeatable analysis of HD chip.

Reference	Presentation	Format
25-NS0001	Vial	1 x 2.0 mL

Reagent to prevent sedimentation of samples that will be tested using the HD chip for T-TAS®01.

Components

- 1 vial x 2.0 mL of Anti-Sedimentation Reagent

