

HUMAN PLASMAS

Normal donor serum

NORMAL SERUM POOL

Fresh frozen serum

Pool of fresh serum from healthy donors



Associated products

Normal donor serum

Informations

The serum is freed from coagulation factors and fibrinogen.

It is obtained by sampling on dry tubes without anticoagulant.

Reference	Presentation	Format
6-SPOOL	Kit	10 x 1.0 mL
6-SPOOL-350	Kit	10 x 0.35 mL

Pool of fresh frozen normal human sera.

The serum pool is collected with great care from healthy male and female donors without drug treatment between 18 and 66 years old. The result is a very high quality product.

Components

- 10 cryotubes x 0.35 mL or 1 mL

Advantages

- Normal human serum, pool of at least 20 sera from at least 20 different healthy donors, decanted, centrifuged and frozen within 3 hours of collection.
- Packaging in plastic cryotubes.

Characteristics

- No additives or preservatives
- No reconstitution error
- Ready to use after thawing (3 min at 37 °C) for 1 mL tubes
- This plasma is stable, if stored at -40 to -80 °C, until the end of the month of the expiration date indicated on the package
- Quality control : example : dosage of the complement

